


1
After threaded implant placement, attach transfers with 1.25mmD Hex Tool. For implants with a Fixture Mount/Transfer already attached, proceed to step 2. Optional: Long impression screws may be used for open-tray impression technique.



2
Block out the top of the transfer and any exposed sutures. Place light body impression material around each transfer and record a full-arch impression with standard body material. Remove the impression after it fully sets. Unscrew transfers from implants and forward with impression for fabrication of a stone model.



3
Connect the transfer to a corresponding replica and insert assembly back into the impression. Create stone model. If desired, use soft tissue material to represent the gingival tissues. Remove transfers.



4
Attach Bar Gold Copings to replicas with a 1.25mmD Hex Tool.
Note: Alternative technique involves use of screw-retained abutments with gold copings for cast bars.



5
Cut and shape gold bar patterns to fit between the Bar Gold Copings. Lute them in place with an autopolymerizing resin parallel to the plane of occlusion. Remove framework from the stone model. Solder the bar segments to the Bar Gold Copings. Finish and lightly polish the soldered bar. Caution: Do not over-polish the areas of the bar that will accept the clip(s).



6
Verify the bar fits passively in the patient's mouth. If necessary, correct the bar by sectioning and resoldering.



7
There are two alternatives for processing the clips to the bar: laboratory or intraoral. For the laboratory technique, remove the bar and send back to the lab for processing of the clips on the stone model. The intraoral technique allows the dentist to pick up the clips off the bar. The patient's existing denture can be used or a new denture fabricated by the laboratory. To pick up the clips, follow traditional prosthodontic techniques.



8
Secure the bar in place with fixation screws and tighten to the proper torque with a calibrated prosthetic torque wrench. Reseat the denture and make final adjustments.

*Immediate loading is indicated for use with Tapered Screw-Vent®, AdVent®, Screw-Vent®, SwissPlus® and Spline Twist® implant systems only. For use with Zimmer Dental AdVent, Tapered Screw-Vent, Spline Twist, Screw-Vent, and SwissPlus implant systems in fully edentulous mandibles.

Tapered Screw-Vent® Implants

Qty.	MTX™			
	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	TSVB8	3.7mm	3.5mm	8mm
<input type="radio"/>	TSVB10	3.7mm	3.5mm	10mm
<input type="radio"/>	TSVB13	3.7mm	3.5mm	13mm
<input type="radio"/>	TSVB16	3.7mm	3.5mm	16mm
<input type="radio"/>	TSVWB8	4.7mm	4.5mm	8mm
<input type="radio"/>	TSVWB10	4.7mm	4.5mm	10mm
<input type="radio"/>	TSVWB13	4.7mm	4.5mm	13mm
<input type="radio"/>	TSVWB16	4.7mm	4.5mm	16mm
<input type="radio"/>	TSV6B8	6.0mm	5.7mm	8mm
<input type="radio"/>	TSV6B10	6.0mm	5.7mm	10mm
<input type="radio"/>	TSV6B13	6.0mm	5.7mm	13mm
<input type="radio"/>	TSV6B16	6.0mm	5.7mm	16mm

Qty. MP-1® HA Dual Transition Selective Surface

	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	TSVH8	3.7mm	3.5mm	8mm
<input type="radio"/>	TSVH10	3.7mm	3.5mm	10mm
<input type="radio"/>	TSVH13	3.7mm	3.5mm	13mm
<input type="radio"/>	TSVH16	3.7mm	3.5mm	16mm
<input type="radio"/>	TSVWH8	4.7mm	4.5mm	8mm
<input type="radio"/>	TSVWH10	4.7mm	4.5mm	10mm
<input type="radio"/>	TSVWH13	4.7mm	4.5mm	13mm
<input type="radio"/>	TSVWH16	4.7mm	4.5mm	16mm
<input type="radio"/>	TSV6H8	6.0mm	5.7mm	8mm
<input type="radio"/>	TSV6H10	6.0mm	5.7mm	10mm
<input type="radio"/>	TSV6H13	6.0mm	5.7mm	13mm
<input type="radio"/>	TSV6H16	6.0mm	5.7mm	16mm

Qty. AdVent® Prosthetics, 4.5mmD Platform

<input type="radio"/>	AVIT/4	Indirect Transfer
<input type="radio"/>	AVR	AdVent Replica
<input type="radio"/>	AVGC3	Bar Gold Coping, 3mmL
<input type="radio"/>	AVGC5	Bar Gold Coping, 5mmL
<input type="radio"/>	AVACT	Tapered Abutment, No Added Height
<input type="radio"/>	AVACT3	Tapered Abutment, 3mmL

Note: 5.7mmD platform prosthetic components are compatible with the AdVent 5.7mmD platform.

© 2004 Zimmer Dental Inc. All rights reserved. 4899, Rev. 12/04.

Prosthetics

Qty.	Replicas	
<input type="radio"/>	IA3	3.5mmD
<input type="radio"/>	IA4	4.5mmD
<input type="radio"/>	IA5	5.7mmD

Qty. Extension Screw for FM/T

<input type="radio"/>	DHTS	Optional Screw for Open Tray Impression with FMT
-----------------------	------	--

Qty. "Cast-To" Gold Abutments, Non-Engaging

<input type="radio"/>	NEA3G	3.5mmD
<input type="radio"/>	NEA4G	4.5mmD

Qty. Tapered Abutments, Non-Engaging

	Platform Diameter	Abutment Length	
<input type="radio"/>	TAC2	3.5mm	2mm
<input type="radio"/>	TAC3	3.5mm	3mm
<input type="radio"/>	TAC4	3.5mm	4mm
<input type="radio"/>	TAC5	3.5mm	5mm
<input type="radio"/>	TACW2	4.5mm	2mm
<input type="radio"/>	TACW3	4.5mm	3mm
<input type="radio"/>	TACW4	4.5mm	4mm
<input type="radio"/>	TACW5	4.5mm	5mm
<input type="radio"/>	TA5C2	5.7mm	2mm
<input type="radio"/>	TA5C3	5.7mm	3mm
<input type="radio"/>	TA5C4	5.7mm	4mm

Qty. Tapered Abutment Gold Copings and Bars

<input type="radio"/>	TGC3	Bar Gold Coping, 3mmL
<input type="radio"/>	TGC5	Bar Gold Coping, 5mmL
<input type="radio"/>	ACTGC	Tapered Abutment Gold Coping
<input type="radio"/>	HGB	1.8mm Round Gold Bar Kit, 50mmL
<input type="radio"/>	DGB	Dolder Resilient Gold Bar, 50mmL

AdVent Implants

Qty.	MTX			
	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	AVB8	3.7mm	4.5mm	8mm
<input type="radio"/>	AVB10	3.7mm	4.5mm	10mm
<input type="radio"/>	AVB13	3.7mm	4.5mm	13mm
<input type="radio"/>	AVB16	3.7mm	4.5mm	16mm
<input type="radio"/>	AVWB8	4.7mm	4.5mm	8mm
<input type="radio"/>	AVWB10	4.7mm	4.5mm	10mm
<input type="radio"/>	AVWB13	4.7mm	4.5mm	13mm
<input type="radio"/>	AVWB16	4.7mm	4.5mm	16mm
<input type="radio"/>	AV6B8	4.7mm	5.7mm	8mm
<input type="radio"/>	AV6B10	4.7mm	5.7mm	10mm
<input type="radio"/>	AV6B13	4.7mm	5.7mm	13mm
<input type="radio"/>	AV6B16	4.7mm	5.7mm	16mm

Qty. MP-1 HA Dual Transition Selective Surface

	Implant Diameter	Platform Diameter	Implant Length	
<input type="radio"/>	AVH8	3.7mm	4.5mm	8mm
<input type="radio"/>	AVH10	3.7mm	4.5mm	10mm
<input type="radio"/>	AVH13	3.7mm	4.5mm	13mm
<input type="radio"/>	AVH16	3.7mm	4.5mm	16mm
<input type="radio"/>	AVWH8	4.7mm	4.5mm	8mm
<input type="radio"/>	AVWH10	4.7mm	4.5mm	10mm
<input type="radio"/>	AVWH13	4.7mm	4.5mm	13mm
<input type="radio"/>	AVWH16	4.7mm	4.5mm	16mm
<input type="radio"/>	AV6H8	4.7mm	5.7mm	8mm
<input type="radio"/>	AV6H10	4.7mm	5.7mm	10mm
<input type="radio"/>	AV6H13	4.7mm	5.7mm	13mm
<input type="radio"/>	AV6H16	4.7mm	5.7mm	16mm

From the U.S.

Tel: 800 854 7019
Fax Orders to: 888 225 2483

From Australia

Tel: +61 2 4297 4182 or 1 800 241 916
Fax Orders to: +61 2 429 74190

From Canada

Tel: +1 800 265 0968 or 1 905 567 2073
Fax Orders to: 1 905 567 2076

From France

Tel: +33 (0)1 45 12 35 35
Fax Orders to: +33 (0)1 45 60 04 88

From Germany

Tel: +49 (0)761 4584 722/723
Fax Orders to: +49 (0)761 4584 729

From Israel

Tel: +972 (0)3 6124242
Fax Orders to: +972 (0)3 6124243

From Spain

Tel: +34 93 846 05 43
Fax Orders to: +34 93 845 43 25

Others

Tel: +1 760 929 4300
Fax Orders to: +1(760) 431 7811

Customer Account _____

Ship To:

Dr. Name _____

Phone (_____) _____

Address _____

City _____

State _____ Zip _____

Order Date _____

Ordered By _____

Ship Via Federal Express

2 Day Overnight (5:30 p.m.) P1 (10:30 a.m.)

Must Arrive By _____

Payment Options

Invoice COD Visa MasterCard AmEx

Credit Card # _____

Exp. Date _____

Signature _____